AT04003, AT04005, AT04006
TEMLAK INSTRUKCJA UŻYTKOWANIA
ARM SLING INSTRUCTION MANUAL
ARMSCHLINGE BEDIENUNGSANLEITUNG
NÁVOD K POUŽITÍ ZÁVĚS RAMENNÍHO KLOUBU PRODYŠNÝ
NÁVOD NA POUŽITIE RAMENNÉHO POPRUHU
MANUEL D'INSTRUCTION POUR L'ÉCHARPE DE BRAS
HANDLEIDING ARM SLING
ESLINGA DE BRAZO MANUAL DE INSTRUCCIONES
MANUALE D'ISTRUZIONE DELL'IMBRAGATURA A BRACCIO
INSTRUKTIONSMANUAL FÖR ARMSTÖD

MD

PL

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ENG

APPLICATION

The arm sling is used to stabilise the upper limb in mild sprains, as a dressing that supports plaster (gypsum), during rehabilitation as well as in case of dislocations of the elbow joint.

USE

Insert the diseased limb into the sling, and put the strap on which it is suspended over the head, placing it on the healthy arm. The length of the strap must be adjusted so that the arm is bent at the elbow joint. AT04006 - in this model you need to fasten an additional Velcro strap.

CONTRAINDICATIONS

physical or mental limitations (e.g. visual impairment) that prevent safe handling.

SIZES

Please measure the length of your forearm to find the correct size:

XXS: (19 - 20.7cm)

XS: (21 - 23 cm)

S: (23.5 - 25.4 cm)

M: (25.7 - 30.5 cm)

L: (30.8 - 35.6 cm)

XL: (35.9 - 40.6 cm)

WASHING INSTRUCTIONS:

The product can only be washed by hand in lukewarm water (up to 40°C) with a mild detergent (e.g. grey soap). Do not iron, do not dry clean, do not use bleach, do not spin dry. Dry unfolded, squeeze gently before drying.

STORAGE

Store in a dry place, away from sunlight.

ATTENTION

In the event of a device-related "serious incident" that has directly or indirectly led to, could have led to or is likely to lead to any of the following:

- (a) death of the patient, user or other person or
- (b) temporary or permanent deterioration of the health of a patient, user or other person, or
- (c) a serious risk to public health

the above 'serious incident' must be notified to the manufacturer and to the competent authority of the Member State in which the user or patient is resident. In the case of Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

NOTES

In case of pain, allergic reactions or other disturbing symptoms unclear to the user related to the use of the medical device, consult a healthcare professional.

NOTE

It is forbidden to use the product other than in accordance with its intended use.

HOW TO DISPOSE OF THE PRODUCT AT END-OF-LIFE

After the device has been withdrawn from service, the medical device can be disposed of as normal municipal waste.

TARGET PATIENT GROUP

People suffering from diseases, dysfunctions or injuries for the treatment/rehabilitation or compensation of which this device is intended (see section on intended use of the device in these instructions). The device can be purchased by the user on their own or on the recommendation of a doctor, therapist or other specialist. Whether you purchase the device yourself or on the recommendation of a doctor/therapist/other specialist, you must take into account the available sizes/necessary functions/sizes/variants of the device, the indications and contraindications for use, as well as the information provided by the manufacturer.

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Model:
Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. AN TAR MEDICAL Sp. Z o.o., 03-068 Warsaw, ul. Zawiślarska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period AN TAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. AN TAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to AN TAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR MEDICAL Sp. z o.o.
03-068 Warsaw, ul. Zawiślarska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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