AT04601 , AT04602, AT04603, AT04604, AT04605

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USER MANUAL FOR BELTS

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NÁVOD NA POUŽITIE BRUŠNÉ PÁSY

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AT04601 Abdominal belt (white)

The abdominal belt is used to strengthen the abdominal muscles and prevent displacement of internal organs caused by abdominal hernias. Velcro fastening gives a wide possibility of adjustment. Recommended in case of abdominal hernia, abdominal muscle weakness, postpartum and after abdominal surgery or surgical procedures.

Belt heights available: 24 cm, 30 cm

Size (waist circumference): XS (<80 cm), S (80 ÷ 90 cm), M (91 ÷ 100 cm), L (101 ÷ 110 cm), XL (111 ÷ 120 cm), XXL (121 ÷ 140 cm), XXXL (> 140 cm)

AT04602 Abdominal belt (black, beige)

The abdominal belt is used to strengthen the abdominal muscles and prevent displacement of internal organs caused by abdominal hernias. Velcro fastening gives a wide possibility of adjustment. Recommended in case of abdominal hernia, abdominal muscle weakness, postpartum and after abdominal surgery or surgical procedures.

Belt Height: 24 cm,

Size (waist circumference): XS (51,6 ÷ 61,6 cm), S (61,6 ÷ 81,3 cm), M (81,3 ÷ 94 cm), L (94 ÷ 104,1 cm), XL (104,1 ÷ 114,3 cm), XXL (114,3 ÷ 121,9 cm), XXXL (122 ÷ 140 cm)

AT04603 Abdominal belt - ostomy

The abdominal-stomal belt is used to support, stabilize and hide the ostomy bag, as well as supporting and strengthening the abdominal wall. It has a zone where the user can independently cut a hole for the ostomy bag - according to your individual needs. Velcro closure gives a wide range of adjustment, and a Velcro-fastened fabric element allows you to hide the ostomy bag. It is made of a soft, breathable material, ensuring great comfort of use. Recommended in case of a stoma, abdominal hernia, abdominal muscle weakness, after abdominal cavity operations, after surgery.

Belt height: 24 cm

Size (waist circumference): XS (70 ÷ 80 cm), S (81 ÷ 90 cm), M (91 ÷ 100 cm), L (101 ÷ 110 cm), XL (111 ÷ 120cm), XXL (121 ÷ 130 cm), XXXL (> 130 cm)

The belt should be positioned so that the ostomy bag opening is positioned over the fistula intestinal. The ostomy bag attached to the belt should be secured with a covering material. Note: The package does not include ostomy bags.

AT04604 Abdominal belt with pad

The abdominal belt is used to strengthen the abdominal muscles and prevent displacement of internal organs caused by abdominal hernias. Velcro fastening gives a wide possibility of adjustment. The attached pad allows you to provide additional pressure in place of the abdominal hernia. The belt is made of soft, breathable material, providing great comfort of use. Recommended in case of abdominal hernia, weakness of the abdominal muscles after the past abdominal surgery, after surgical procedures.

Belt height: 24 cm

Size (waist circumference): S (61,6 ÷ 81,3 cm), M (81,3 ÷ 94 cm), L (94 ÷ 104,1 cm), XL (104,1 ÷ 114,3 cm), XXL (114,3 ÷ 121,9 cm), XXXL 121,9 – 132,9 cm

AT04605 Abdominal belt

The abdominal belt is used to strengthen the abdominal muscles and prevent displacement of organs internal due to abdominal hernias. It is made of a unique composition of two meshes, outer - strong, carrying loads and inner - soft and comfortable to touch with the body. Thanks to this, the belt is extremely breathable, provides excellent drainage of excess heat and sweat, and the „breathing“ of the body beneath it. Velcro fastening gives a wide possibility of adjustment. Recommended in the case of abdominal hernia, abdominal muscle weakness, postpartum and after abdominal surgery or surgical procedures.

Belt height: 24 cm

Sizes: S (70 ÷ 90 cm), M (80 ÷ 95 cm), L (90 ÷ 105 cm), XL (95 ÷ 112 cm), XXL (105 ÷ 121 cm), XXXL (115 ÷ 140 cm)

WARNING:

In case of a device-related "serious incident", which directly or indirectly led to, could have led to, or is likely to lead to any of the following events:

- death of a patient, user or another person, or
- temporary or permanent deterioration of the health of a patient, user or another person, or
- a serious risk to public health, such a "serious incident" must be notified to the manufacturer and to the competent authority of the Member State, in which the user or patient resides. For Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

WARNING:

In the event of pain, allergic reactions or other alarming, unclear symptoms related to the use of the medical device, please consult a healthcare professional.

User manual:

In order to select the correct size, measure the circumference of the abdomen.

Undo the Velcro. Put the belt on the lumbosacral spine, and then fasten it with Velcro at the front.

The product has a wide range of adjustments. The product can be quickly and easily adjusted to the shape of the patient's body. Velcro fastening guarantees effective immobilization.


Fig. 1

Fig. 2

Fig. 3

Fig. 4

Comments and recommendations:

Please read this manual carefully before putting the product on.

Before each use, make sure that the product does not have any damages at the seams and other elements.

Using the product contrary to the rules contained in the instructions may result in serious complications.

The product should be used according to the recommendations of a doctor or physiotherapist.

The product should not be used in the case of open wounds, abrasions of the epidermis, as well as after applying a warming gel or ointment.

Wearing the product during the night is only permitted after the recommendation of your treating physician.

The condition of the skin in the place where the product adheres to, should be checked. If there are skin changes contact your doctor.

· Unless your doctor recommends otherwise, you should wear the device no longer than 4 hours a day.

NOTE: It is forbidden to use the product contrary to its intended use!

TARGET PATIENT GROUP

Persons suffering from diseases, dysfunctions or injuries for the treatment/rehabilitation or compensation of which this device is intended (see section on intended use of the device in these instructions). The device can be purchased by the user on their own or on the recommendation of a doctor, therapist or other specialist. Whether you purchase the device yourself or on the recommendation of a doctor/therapist/other specialist, you must take into account the available sizes/necessary functions/size/variants of the device, the indications and contraindications for use, as well as the information provided by the manufacturer.

CONTRAINDICATIONS

Physical or mental limitations (e.g. visual impairment) that prevent the safe handling of the product.

CLEANING AND MAINTENANCE

The brace can only be washed by hand in lukewarm water (at the temperature up to 40 °C) with a mild detergent (e.g. gray soap).

Do not iron. Do not dry clean. Do not use bleach. Protect against moisture. Do not expose to high temperatures. Do not tumble dry. Dry open product, squeeze gently before drying.

AT04604-**WARNING:** The product contains buckwheat hulls. The orthopedic pad cannot be washed.

HOW TO DISPOSE OF THE PRODUCT

The medical device can be disposed of as normal municipal waste.

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Größe (Taillenumfang): XS (70 ÷ 80 cm), S (81 ÷ 90 cm), M (91 ÷ 100 cm), L (101 ÷ 110 cm), XL (111 ÷ 120 cm), XXL (121 ÷ 130 cm), XXXL (> 130 cm)

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Größen: S (70 ÷ 90 cm), M (80 ÷ 95 cm), L (90 ÷ 105 cm), XL (95 ÷ 112 cm), XXL (105 ÷ 121 cm), XXXL (115 ÷ 140 cm)

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WARRANTY CARD

Model:
Serial number:

Point of sale stamp
and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawisłańska 43, hereby grants a 12-month warranty on the rollator, starting from the date of purchase by the Buyer. Products financed by the National Health Fund are covered with the extended warranty period, which covers half of the useful life specified in the Regulation on Guaranteed Services, provided that the useful life is specified in the Regulation in months or years.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, send it immediately to ANTAR (as specified in point 1) or contact the point of sale.
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum permissible weight of the user as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts subject to normal wear and tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED. IN ORDER TO EXTEND WARRANTY PERIOD FOR PRODUCTS FINANCED BY THE NATIONAL HEALTH FUND, ATTACH A COPY OF THE ORDER FOR THE SUPPLY OF MEDICAL DEVICES BEING ORTHOPAEDIC ITEMS AND AIDS.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

 ANTAR Sp. J.
03-068 Warsaw, ul. Zawisłańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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