AT51028
CHODZIK INSTRUKCJA UŻYTKOWANIA
USER MANUAL ROLLATOR
BENUTZERHANDBUCH ROLLATOR
UŽIVATELSKÝ MANUÁL CHODÍTKO
NÁVOD NA POUŽITIE CHODÍTKO
MANUEL DE L'UTILISATEUR ROLLATEUR
HANDLEIDING ROLLATOR
MANUAL DE USUARIO ANDADOR
MANUALE D'USO ROLLATORE
ANVÄNDARMANUAL ROLLATOR

PL

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EN

Thank You for Purchasing Our Product

Before using the walker, please read this user manual carefully. Incorrect use may result in injuries, bodily harm, or accidents. To ensure safe and successful use, we strongly recommend thoroughly reviewing this manual.

NOTICE:

Check all parts for any damage that may have occurred during transport. If you notice any damage, DO NOT use the product. For further information, contact the manufacturer.

GENERAL WARNING

DO NOT assemble or use the product before reading and understanding this manual. If you do not understand the warnings, precautions, or instructions, consult a healthcare professional, distributor, or technical specialist before attempting assembly to prevent injury or damage to the product. Every user should consult their physician to determine the appropriate adjustment and usage.

The walker is not intended for the user to sit on and propel themselves forward. A physician or specialist should be present when adjusting the walker's height to provide support and ensure proper tightening of screws. When the walker is fully unfolded, always make sure the tube under the seat is locked.

Ensure that the height adjustment knobs are locked in place so that the main supporting structure does not move under pressure. Failure to follow these instructions may result in injuries. Before use, verify that ALL handle and seat adjustment components are properly tightened and that the wheels and moving parts are in good condition.

During use, ALL wheels must REMAIN IN CONTACT with the ground at ALL TIMES to ensure the walker's stability. Failure to comply may result in the walker tipping over and causing injury to the user.

NOTICE:

In the event of a **serious incident** related to the product that has directly or indirectly resulted, may have resulted, or could result in any of the following:

- a) Death of the patient, user, or another person
- b) Temporary or permanent deterioration of the health of the patient, user, or another person
- c) A serious public health risk

This **serious incident** must be reported to the manufacturer and the competent authority of the member state where the user or patient resides. In Poland, the responsible authority is the **Office for Registration of Medicinal Products, Medical Devices, and Biocidal Products**.

NOTICE:

During use, handling, folding, and adjusting mechanisms, there is a risk of entrapment and/or pinching of body parts between components. These actions should be performed with extreme caution. After adjustments, stabilize the position by firmly tightening nuts/screws.

NOTICE:

If you experience pain, allergic reactions, or any other concerning or unclear symptoms related to the use of this medical device, consult a healthcare professional.

NOTICE:

The manufacturer is not responsible for damages caused by lack of maintenance, improper servicing, or failure to comply with the instructions in this manual.

NOTICE:

The product must **not** be used for any purpose other than its intended use.

INTENDED USER GROUP

Individuals experiencing medical conditions, dysfunctions, or injuries requiring treatment, rehabilitation, or compensation for which this product is designed (see the section on product purpose in this manual). The product may be purchased independently or based on a physician's, therapist's, or other specialist's recommendation. Whether purchasing independently or upon medical advice, users should consider available sizes, necessary functions, product variants, indications, contraindications, and manufacturer-provided information.

PURPOSE

This product is designed to assist individuals with limited mobility. The walker is intended for **single-person** use only.

INDICATIONS

Inability to move independently due to age, illness, or general health condition

Balance disorders

Post-operative recovery

Rehabilitation

CONTRAINDICATIONS

Physical or mental impairments (e.g., vision impairment) that prevent the safe use of the walker.

WARNINGS REGARDING ASSEMBLY AND INSTALLATION

After assembling the walker, ensure that:

The seat is locked in place.

Both handle tubes are securely attached to the main support frame using nuts.

Both handles are at the same height and do not shift.

By following these guidelines, you can ensure safe and effective use of the walker.

Ensure Safe Locking of the Seat Tube

Before use, always check that the locking components are in their proper positions.

UNFOLDING

Set up the walker and lock the tube under the seat.

Attach the basket to the hooks and lower the seat. (Fig. 1)

Insert the handle tubes into the corresponding openings, adjust them to the proper height, and tighten the knobs. (Fig. 2)

Finally, install the backrest, ensuring the copper pins fit into the correct holes. (Fig. 3)

FOLDING

Remove the backrest and handles. Lift the seat, remove the basket, then unlock the tube under the seat to fold the product.

BRAKING

The walker brakes when pressure is applied to the handles.

TECHNICAL SPECIFICATIONS

Type: Four-wheeled, foldable aluminum walker

Wheel size: 5"

Total height: 78 cm – 88 cm

Total length: 69 cm

Upper width: 50 cm

Total width: 58 cm

Seat thickness: 4 cm

Seat width: 36 cm

Seat height from ground: 54 cm

Seat depth: 31 cm

Weight: 5.75 kg

Maximum user weight: 120 kg

Distance between handles: 45 cm

THIS SYMBOL INDICATES THE MAXIMUM USER WEIGHT.

PACKAGE CONTENTS

Walker

Basket

Two black adjustment knobs

Two handle tubes

Backrest

User manual

MAINTENANCE / CLEANING

Ensure that all components are always properly attached.

Check all fastening screws and tighten them if necessary.

Immediately replace any damaged or worn-out parts.

Clean the product with a damp cloth and mild detergent, then wipe it dry.

Do not use abrasive cleaning tools.

WHEEL MAINTENANCE

Clean the wheels with warm water and a mild cleaning agent. Dry them with a cloth.

Periodically lubricate with **3-in-1 oil** if the wheels begin to squeak.

STORAGE

Store the product in a dry, cool place, away from moisture and direct sunlight.

DISPOSAL OF THE MEDICAL DEVICE

Dispose of the medical device in accordance with applicable disposal regulations.

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Model:
Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. AN TAR MEDICAL Sp. z o.o., 03-068 Warsaw, ul. Zawisłańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period AN TAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. AN TAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to AN TAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR MEDICAL Sp. z o.o.
03-068 Warsaw, ul. Zawisłańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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