AT51034

**INSTRUKCJA UŻYTKOWANIA CHODZIK
USER MANUAL ROLLATOR
GEBRAUCHSANWEISUNG ROLLATOR
NÁVOD K POUŽITÍ CHODÍTKO ČTYŘKOLOVÉ
NÁVOD NA POUŽITIE CHODÍTKO
MANUEL D'UTILISATION ROLLATOR
GEBUIKERSHANDLEIDING ROLLATOR
MANUAL DEL USUARIO ROLLATOR
MANUALE D'USO ROLLATOR
ANVÄNDARHANDBOK ROLLATOR**

IU_AT51034

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NOTE: Inspect all parts for shipping damage. If any such a damage is found, DO NOT use the product. For more information please contact the manufacturer.

RECOMMENDATIONS SAFETY

To be safe when using the rollator, please follow these guidelines:

GENERAL WARNING

Do NOT assemble and use the product until you have read and understood this manual. If you do not understand any warnings, cautions, or instructions, contact your healthcare professional, distributor, or service technician before attempting to assemble the product to avoid personal injury or damage to the product. Each user should ask his physician or therapist to determine the appropriate variant of regulation and method of use. The rollator is not intended to be guided by the user while sitting on it. A physician / therapist should be present when adjusting the height of a rollator to provide support and support for the user

setting the brakes and tightening the screws. Always make sure the seat tube is locked and properly attached to the frame when the rollator is unfolded, and that the brake cable is located inside the side tube and is tightly attached to it. Check after each height adjustment the brake cable to avoid loosening the cable. Make sure the knobs are locked so that the main support structure it will not move under pressure. The quick release mechanism is installed in place and will remain in the release position even if you pull forcefully in the opposite direction. The copper pins in the regulating tube pop out on the right. Failure to follow the above recommendations may lead to injuries. Before use, make sure that ALL the adjustment elements of the handles and the seat are secure and

that the wheels and moving parts of the product are in good condition. During use, ALL the wheels must be in contact with the ground AT ALL TIMES to ensure the stability of the rollator. BEFORE the user is seated on the seat, the brakes MUST be in the locked position.

When the rollator is used in a stationary position, its handbrakes MUST be in the locked position.

WARNING:

In case a serious incident that has occurred in relation to the medical device, that directly or indirectly led, might have led or might lead to any of the following:

- the death of patient, user or other person,
 - the temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
 - a serious public health threat;
- the above mentioned serious incident should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is resident;

WARNING:

In case of any pain, allergic reaction or other worrying symptoms, unclear for the user which have occurred in relation to the medical device, the user

should consult a healthcare professional;

NOTE: The manufacturer is not responsible for damage caused by neglect of improper maintenance, servicing or resulting from non-compliance with the recommendations contained in this manual.

ATTENTION: it is forbidden to use the product other than in accordance with its intended purpose!

NOTE: There may be a risk of tipping over if the product is used improperly. Please follow the sitting / standing / moving instructions. After completing the adjustment, stabilize the position by tightening the nuts / bolts carefully.

WARNINGS PRODUCT STABILITY

The maximum load for the rollator is 136 kg, INCLUDING the contents of the bag. The maximum load of the rollator bag is 5 kg.

WARNINGS ASSEMBLY AND ASSEMBLY

After disassembling or assembling the rollator, make sure the seat is locked in both tubes and both grip tubes are locked in the main support frame with bolts and nuts. Also, make sure that the copper pins come out of the corresponding holes in the tubes of the wheels and that they are properly locked, and that both handles are at the same height and do not move.

WARNINGS INSTALLATION

Make sure the seat tube is securely fastened to the lock, all copper pins slide out of their respective holes, and the front and rear tubes are securely inserted into the main frame and secured with knobs. BEFORE use, always check that the rollator and its components are securely locked in place.

INTENDED USE: The product is intended to support the movement of people with reduced walking ability. The rollator is designed for one person.

CONTRAINDICATIONS: physical or mental limitations (e.g. visual impairment) that prevent safe handling with the rollator.

DELIVERED PARTS

Rollator, bag, instruction manual

UNFOLDING

Place the rollator on an even surface, attach the wheels, then unfold the rollator and mount the bag.

Fig. 1

Fig. 2

Fig. 3

ADJUSTING TUBE HANDLE

Press the latch on the adjustment tube and slide it up or down to adjust the height. Make sure you adjust the latch after adjusting it in the appropriate hole.

FOLDING

Lift the handle of the seat to fold the rollator as shown in fig. 4

Fig. 4


Fig. 5

APPLICATION / LOCKING / UNLOCKING OF HAND BRAKES

NOTE: For more information, see WARNINGS. INSTALLATION IN THE RECOMMENDATIONS SAFETY OF THIS MANUAL.

CAUTION: When using and handling the rollator, and when folding it and adjusting the mechanisms, there may be a risk of entrapment and / or pinching of body parts of the user / accompanying person in openings / in gaps between the components. These steps must be performed with extreme caution. After completing the adjustment, stabilize the position by accurately tightening nuts / bolts.

APPLICATION OF HAND BRAKES

1. Pull the brake towards the handgrips.
2. Do one of the following:
 - A. STANDING - hold the brake handle up
 - B. MOVING - Release the brake handle

LOCKING THE HAND BRAKES. UNLOCKING THE HAND BRAKES

Push down on the bottom of the brake handle as shown until you hear a click (Fig. 8).

UNLOCKING THE HAND BRAKES

Pull the brake handles up until it is released as shown in fig. 5. Release the brake handle.

Make sure the brake cable runs along the sideline and the unused part is coiled in front of the rollator.

Fig. 7

Fig. 8

MAINTENANCE

1. Please ensure that all attachments are always secured.

2. Check the bolt securing the seat to the side elements. Tighten as needed.

3. Replace any damaged or damaged parts immediately.

WHEEL MAINTENANCE

1. Wheels should be cleaned with warm water and a mild detergent. Dry with a dry cloth.
2. Periodically lubricate with 3-in-1 grease if the wheels begin to squeak.

TECHNICAL PARAMETERS

- maximum height: 98 cm
- minimum height: 87 cm
- width between the handles: 50 cm
- seat height: 55 cm
- seat width: 45 cm
- seat depth: 24 cm
- maximum width: 56 cm
- overall length: 85 cm
- width after folding: 27 cm
- wheel diameter: 19 cm
- wheel thickness: 4 cm
- 2 shock absorbers
- maximum load of the bag: 5 kg
- maximum load: 136 kg

This symbol indicates the user's maximum weight

SAFE DISPOSAL OF THE MEDICAL DEVICE

When the device has been withdrawn from use it can be treated as municipal waste except for electric medical devices – in such case the electric medical devices shall be disposed in the way described by the law for electrical and electronic waste disposal.

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WARRANTY CARD

Model:
Serial number:

Point of sale stamp
and a legible signature of the Seller:

1. AN TAR Sp. J., 03-068 Warsaw, ul. Zawisłarska 43, hereby grants a 12-month warranty on the rollator, starting from the date of purchase by the Buyer. Products financed by the National Health Fund are covered with the extended warranty period, which covers half of the useful life specified in the Regulation on Guaranteed Services, provided that the useful life is specified in the Regulation in months or years.
2. During the warranty period AN TAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. AN TAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, send it immediately to AN TAR (as specified in point 1) or contact the point of sale.
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum permissible weight of the user as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts subject to normal wear and tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED. IN ORDER TO EXTEND WARRANTY PERIOD FOR PRODUCTS FINANCED BY THE NATIONAL HEALTH FUND, ATTACH A COPY OF THE ORDER FOR THE SUPPLY OF MEDICAL DEVICES BEING ORTHOPAEDIC ITEMS AND AIDS.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawisłarska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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