## **AT51035**

**CHODZIK INSTRUKCJA UŻYTKOWANIA  
ROLLATOR USER MANUAL  
ROLLATOR-BENUTZERHANDBUCH  
UŽIVATELSKÁ PŘÍRUČKA K ROLÁTORU  
POUŽÍVATELSKÁ PRÍRUČKA K ROLÁTORU  
MANUEL D'UTILISATION DU ROLLATOR  
GEBRUIKSAANWIJZING ROLLATOR  
ROLLATOR MANUAL DE USUARIO  
MANUALE D'USO DEL ROLLATOR  
ROLLATOR ANVÄNDARMANUAL**

**MD**

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**ENG**

#### **RECOMMENDATIONS SAFETY**

To stay safe while using the walker, you **MUST** follow these guidelines:

**NOTE:** Inspect all parts for shipping damage. If you notice such damage, **DO NOT** use the product. More information from the manufacturer.

#### **GENERAL WARNING**

**DO NOT** assemble and use the product before reading and understanding this manual. If you do not understand any warnings, cautions, or instructions, contact your healthcare professional, distributor, or technician before attempting to assemble the product to avoid injury or damage to the product. Each user should ask his doctor or therapist to determine the appropriate adjustment variant and method of use. The walker is not intended to be controlled by the user sitting on it. A physician/therapist should be present when adjusting the height of the walker to ensure the user is supported and that the brakes are set correctly and the screws are tightened. When the walker is unfolded, always make sure that the tube at the seat is locked and properly attached to the frame and that the brake cable is inside the side tube and is firmly attached to it. After each height adjustment, check the brake cable to avoid cable slack. Make sure that the knobs are locked in such a way that the main supporting structure will not move due to pressure. The copper pins in the regulating tube extend outwards. Failure to follow the above instructions may result in injury. Before use, check that **ALL** handle and seat adjustment elements have been secured and that the wheels and movable parts of the product are in good condition. During use, **ALL** the wheels must **STAY** in contact with the ground at **ALL TIMES** to ensure the stability of the walker. **BEFORE** the user sits on the seat, the brakes **MUST** be in the lock position. When the walker is used in a stationary position, the hand brakes **MUST** be in the lock position.

#### **ATTENTION:**

In the event of pain, allergic reactions or other disturbing, unclear for the user symptoms related to the use of the medical device, consult a healthcare professional.

#### **ATTENTION:**

In the event of a product-related "serious incident" that directly or indirectly resulted in, lead to or may lead to any of the following events:

- a) death of a patient, user or other person or
- b) temporary or permanent deterioration of the patient's, user's or other person's health or
- c) a serious threat to public health

this "serious incident" should be reported to the manufacturer and to the competent authority of the Member State where the user or patient resides. In the case of Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

**NOTE:** The manufacturer is not responsible for damage caused by neglected maintenance, improper servicing or as a result of failure to comply with the cleaning rules or recommendations contained in this manual.

**ATTENTION:** It is forbidden to use the product in a way other than for its intended purpose!

**NOTE:** If the product is used incorrectly, there may be a risk of tipping over. Please follow the getting on/getting off/moving instructions. After adjustment, stabilize the position by carefully tightening the nuts/bolts.

#### **SAFETY WARNINGS PRODUCT STABILITY**

The maximum load on the walker is 136 kg, **INCLUDING** the contents of the bag. The maximum load of the walker bag is 5 kg.

#### **SAFETY WARNINGS ASSEMBLY AND DISASSEMBLY**

After disassembling or assembling the walker, make sure that the seat is locked in both tubes, and both handle tubes - in the main support frame, with bolts and nuts. Also make sure that the copper pins come out of the corresponding holes in the caster tubes and are properly locked, and that both handles are at the same height and do not move.

#### **SAFETY WARNINGS INSTALLATION**


Make sure the seat tube is securely fastened in the lock, all copper pins slide out of their respective holes, and the front and rear tubes are securely inserted into the main frame and secured with knobs. Always check that the walker and its components are securely locked in place **BEFORE** use.

**DESTINATION:** The product is intended to support people with limited walking ability while moving. The walker is designed for one person.

**CONTRAINDICATIONS:** physical or mental limitations (e.g. visual impairment) that prevent safe handling of the walker.

### UNFOLDING

1. Place the walker on a level surface, unfold the main frame and then place the handles in the frame so that the copper pins are in the corresponding holes outside. The brake cable should always be on the inside of the frame. (fig.1)
2. Attach the shock absorbers to the main frame using the round hooks. (fig.2)
3. Attach the rear wheels to the main frame using the square hooks. (fig.3)
4. Secure the front wheels with the black bolts, pads and nuts. (fig.4 and 5)
5. Place the bag on the front of the walker and the cane holder in place. (fig.6 and 7)

### HANDLE TUBE POSITIONING

Press the latch (pin) in the adjustment tube and move it up or down to adjust the height. Make sure that after adjustment the latch (pin) is in the correct hole and the copper pin is extended outwards. (fig.1)

### FOLDING

To completely disassemble the walker, perform all the above steps in reverse order.

You can also fold the walker only for transport, just pull the handle on the seat and the walker will be folded.

### APPLYING/LOCKING/UNLOCKING THE HANDBRAKES


ATTENTION: During the use and handling of the walker, as well as when folding it and adjusting the mechanisms, there may be a risk of entrapment and/or pinching of the user's/assistant's body parts in the openings/gaps between the elements. These steps should be performed with particular care. After adjustment, stabilize the position by carefully tightening the nuts/bolts.

### USING THE HANDBRAKES

1. Pull the brake towards the hand grips.
2. Do one of the following:
  - A. WITHOUT MOVING - keep the brake handle up:
  - B. MOVING - release the brake handle

### LOCKING THE HANDBRAKES

1. Push down the bottom of the brake handle as shown until you hear a click (Fig.7)

6 (unlocked state) 7 (locked state)

### RELEASING THE HANDBRAKES

Pull the brake handles up until they release as shown in Fig. 6.

1. Release the brake handle.
2. Make sure the brake cable runs along the side tube and the unused part is rolled up at the front of the walker.

**DESTINATION:** The product is intended to support people with limited walking ability while moving.

The walker is designed for one person.

**CONTRAINDICATIONS:** physical or mental limitations (e.g. visual impairment) that prevent safe handling with a walker.

**CHARACTERISTICS:**

- Lightweight, aluminum
- Possibility of quick unfolding/folding of the walker
- Has a seat function

**TECHNICAL PARAMETERS**

Four-wheel aluminum walker, foldable.

Weight: 8.6 kg

Height: 87-98 cm

Total length: 78 cm

Folding width: 27 cm

Wheel diameter: front 8", rear: 8"

Seat height: 53.5 cm

Width between handles: 50 cm

Seat width: 45 cm

Seat depth: 24 cm

Distance between the front wheels: 50 cm

Distance between front wheel and back wheel: 51 cm


Maximum user weight: 136 kg

**THIS MARK INDICATES THE MAXIMUM USER WEIGHT**

**INCLUDED**

**Walker:** Main frame, front wheels, rear wheels with handles and brakes, bag, cane holder, 2 screws, 2 pads, 2 nuts, 2 round hooks, 2 square hooks.

**TARGET GROUP OF PATIENTS**

The healthcare professional indicates on his own responsibility the use of the device for adults and children, taking into account the available variants/sizes/necessary functions/sizes and indications, bearing in mind the information provided by the manufacturer.

**HOW TO DISPOSAL OF THE PRODUCT**

At the end of life, the medical device can be disposed of as normal household waste

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Model: .....  
Serial number: .....

Point of sale stamp

and a legible signature of the Seller: .....

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawiślańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.  
03-068 Warsaw, ul. Zawiślańska 43  
Tel. 22 518 36 00, Fax 22 518 36 30  
www.antar.net

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Antar Medizin GmbH