AT51036

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USER'S MANUAL ROLLATOR

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NÁVOD NA POUŽITIE

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HANDLEIDING ROLLATOR

MANUAL DE USO ROLLATOR

MANUALE D'USO ROLLATOR

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GENERAL WARNING

DO NOT assemble or use the product before reading and understanding these instructions. If you do not understand the warnings, cautions or instructions, contact your healthcare professional, distributor or technician before attempting to assemble the product to avoid injury or damage to the product. Each user should ask their doctor or therapist to determine the appropriate adjustment option and method of use. The rollator is not intended to be steered by the user while sitting on it.

When adjusting the height of the rollator, the doctor/therapist should be present to ensure that the user is supported and that the brakes are correctly adjusted and the bolts tightened. When the rollator is extended, always make sure that the seat tube is locked and properly connected to the frame and that the brake cable is inside the side tube and firmly attached to it. Check the brake cable after each height adjustment to avoid the cable coming loose. Make sure that the knobs are locked so that the main support structure will not move as a result of pressure. The quick release mechanism is fitted in the correct position and will remain in the release position even if pulled hard in the opposite direction. The copper pins in the adjustment tube extend from the right side. Failure to comply with the above may lead to injury. Before use, check that ALL handle and seat adjustments have been secured and that the castors and moving parts of the product are in good condition. During use ALL castors must be in contact with the ground AT ALL TIMES to ensure the stability of the rollator. BEFORE the user sits on the seat, the brakes MUST be in the locked position.

When the rollator is used in a stationary position, its hand brakes MUST be in the locked position.

NOTE: Inspect all parts for damage that has occurred in transit. If such damage is observed, DO NOT use the product. Contact the manufacturer for more information.

WARNINGS REGARDING FOLDING AND ASSEMBLY

After unfolding or assembling the rollator, make sure the seat is locked into both tubes and both handle tubes are locked into the main support frame with bolts and nuts. Also make sure that the copper pins slide out of the corresponding holes in the caster tubes and that they are properly locked, and that both handles are at the same height and do not move.

UNFOLDING

Ensure that the seat tube is securely fastened in the locking mechanism, that all copper pins slide out of their respective holes, and that the front and rear tubes are securely inserted into the main frame and secured with knobs. BEFORE use, always check that the rollator and its components are securely locked in place. Then, screw in the screws to position the front wheels and also both tubes for the rear wheels.

FOLDING

Lift the seat handle to fold the rollator as shown in fig.5

Fig. 1

Fig. 2

Fig. 3


Fig. 4


Fig. 5

ADJUSTING THE HANDLE TUBE

Press the latch in the adjustment tube and move it up or down adjusting the height. Make sure the latch is in the correct hole after adjustment.

NOTE: When using and operating the rollator and when folding and adjusting the mechanisms, there may be a risk of entrapment and/or compression of body parts of the user/assisted person in the openings/gaps between the components. Take one handle and press down on the seat until all four pegs are unfolded and the seat is flat. Perform these operations with particular care. After adjustment, stabilise the position by carefully tightening the nuts/bolts.

WARNING:

In the event of a device-related "serious incident", which directly or indirectly led to, could have led to, or is likely to lead to any of the following events:

- a) death of a patient, user or another person, or
- b) temporary or permanent deterioration of the health of a patient, user or another person, or
- c) a serious risk to public health, such a "serious incident" must be notified to the manufacturer and to the competent authority of the Member State, in which the user or patient resides. For Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

WARNING:

In the event of pain, allergic reactions or other alarming, unclear symptoms related to the use of the medical device, please consult a healthcare professional.

NOTE: It is forbidden to use the product contrary to its intended use!

UNLOCKING THE HANDBRAKE

1. Pull the brake handles upwards until they release as shown in Fig. 6.
2. Release the brake handle.
3. Make sure the brake cable runs along the side tube and the unused portion is coiled at the front of the rollator.

Fig. 6 (unblocked state)

Fig. 7 (blocked state)

APPLICATION OF HANDBRAKE

1. pull the brake towards the hand grips.
2. perform one of the following actions:
A. WITHOUT MOVING - hold the brake handle up;
B. MOVING - release the brake handle

LOCKING THE HANDBRAKE UNLOCKING THE HANDBRAKE

1. push down on the bottom of the brake handle as shown until you hear a click (Fig. 7)
2. For more information refer to the UNLOCKING HAND BRAKE section of this manual.

MAINTENANCE

1. Ensure all attachments are secure at all times.
2. Check the bolt securing the seat to the side members. Tighten if necessary.
3. promptly replace any damaged or deteriorated components.

WHEEL MAINTENANCE

1. Clean castors with warm water and a mild detergent. Dry with a dry cloth.
2. lubricate periodically with 3-in-1 grease if casters start to squeak.

PURPOSE

Auxiliary and rehabilitation walker supporting upright position and movement of disabled people.

TARGET GROUP OF PATIENTS

Persons suffering from diseases, dysfunctions or injuries for the treatment/rehabilitation or compensation of which this device is intended (see section on intended use of the device in this instruction manual). The device can be purchased by the user on their own or on the recommendation of a doctor, therapist or other specialist. Whether you purchase the device yourself or on the recommendation of a doctor/therapist/other specialist, you must take into account the available sizes/necessary functions and variants of the device, the indications and contraindications for use and the information provided by the manufacturer.

CONTRAINDICATIONS

Physical or mental limitations (e.g. visual impairment) that prevent the safe handling of the product.

Technical characteristics:

- four-wheel walker, aluminium, foldable
- maximum height: 103 cm
- minimum height: 89 cm
- 6 steps of adjustment every 2.5 cm
- width between rear wheels: 55 cm
- width between front wheels: 56 cm
- width of seat: 45 cm
- depth of the seat: 25 cm
- height of the seat: 59 cm
- diameter of wheels: 29 cm, inner tube wheels
- distance between handles: 46 cm
- product weight: 9 kg
- bag included
- space for a crutch/blade
- maximum load: 136 kg
- maximum load of the bag: 5 kg

THIS SIGN INDICATES THE MAXIMUM WEIGHT OF THE USER

DISPOSAL OF THE PRODUCT

The product must be disposed of in accordance with the applicable disposal conditions.

CAUTION: There may be a risk of tipping over if the product is used improperly. Please follow the boarding / dismounting / moving instructions. After completing the adjustment, stabilize the positions by tightening the nuts / bolts carefully.

ATTENTION: The manufacturer is not responsible for damage caused by neglect of improper maintenance servicing or resulting from non-compliance with the recommendations contained in this manual.

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WARRANTY CARD

Model:
Serial number:

Point of sale stamp
and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawiańska 43, hereby grants a 12-month warranty on the rollator, starting from the date of purchase by the Buyer. Products financed by the National Health Fund are covered with the extended warranty period, which covers half of the useful life specified in the Regulation on Guaranteed Services, provided that the useful life is specified in the Regulation in months or years.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, send it immediately to ANTAR (as specified in point 1) or contact the point of sale.
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum permissible weight of the user as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts subject to normal wear and tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED. IN ORDER TO EXTEND WARRANTY PERIOD FOR PRODUCTS FINANCED BY THE NATIONAL HEALTH FUND, ATTACH A COPY OF THE ORDER FOR THE SUPPLY OF MEDICAL DEVICES BEING ORTHOPAEDIC ITEMS AND AIDS.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

 ANTAR Sp. J.
03-068 Warsaw, ul. Zawiańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net e-mail: antar@antar.net

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