
# PNEUMATIK ANTI DECUBITUS AIR CELL CUSHION

**Models: AT52108, AT52109, AT52110**

### Pneumatic anti-decubitus cushion

The pneumatic anti-decubitus cushion allows the user a comfortable and stable sitting position. The Single Room Pneumatic Pressure Relief Pillow is designed to evenly distribute the patient's weight and adjust the pillow to the user's current body position. The pillow is soft, lightweight, easy to carry and clean.

The product allows the sacrum to be fully immersed in the pillow, reducing pressure on the coccyx and sciatica. The cushion prevents pressure points, prevents pain during movement, reduces the risk of damage to skin tissue and optimizes tissue relief.

- Combined air chambers are easy to inflate and adjust.
- A specially designed structure of air-filled spaces allows air to flow between them, optimizing the distribution of body pressure on the pad.
- Each pillow contains a valve that is used to adjust the amount of air in the pillow and therefore the softness of the pillow.

### Indications:

- Patients with hip ulcers.
- Patients with pelvic rotation or pelvic changes.
- Patients with stents.
- People in wheelchairs.

Maximum user weight AT52109: 120 kg. Maximum user weight AT52108, AT52110: 136 kg

### WARNING:

1. Read the instruction manual carefully before using the pillow.
2. Be careful of what the bottom part of the pillow contacts as there may be some shifting phenomenon. Be careful.
3. Do not use sharp tools when touching the cushion.
4. Inflate the cushion with the correct amount of air before use.
5. Keep the pillow away from fire and moisture.
6. Do not bend the pillow when filling air.
7. If any part of the body is left on the cushion for too long, air will escape from that area, reducing the therapeutic effect.
8. Check the cushion daily (pressure regulation)
9. The product should be stored in a dry environment.

### CAPACITY

1. The cushion provides comfort and stability by evenly distributing the pressure from the hips so that the user sinks into it when sitting on this cushion.
2. The cushion is intended for contact use, particularly for patients at high risk of bedsores.
3. Filled gaps make your posture more balanced and stable.
4. The cushion is suitable for people with disabilities, the elderly, office workers, drivers, people who sit for a long time, etc.

### PRODUCT SPECIFICATION:



Product code	Dimensions of the pillow Width Length Height	Product features
AT52108	44 x 38 x 5 cm	- small air spaces with no room for growth - 9 rows with 8 air spaces - Dimensions of the air wings: width 4 cm, length 4 cm - Net Weight: 0.33 kg - Material: PVC
AT52109	45 x 39 x 4 cm	- Airspaces at low altitude - Greater stability due to larger and lower spaces - 7 rows with 6 air spaces - Dimensions of the air wings: width 6 cm, length 6 cm - Net weight: 0.255 kg - Material: PVC
AT52110	45 x 39 x 6 cm	- relatively high airspaces - 7 rows with 6 air spaces - More pleasant feel due to larger and higher air spaces - Dimensions of the air wings: width 6 cm, length 6 cm - Net weight: 0.315 kg - Material: TPU

THIS SIGN INDICATES THE MAXIMUM WEIGHT OF THE USER

Air Wings: Made of PVC and TPU.


1. Ventilation: Each pillow contains a valve that can be used to adjust the amount of air in the pillow and thus the softness of the pillow.
2. Instruction Manual: Contains instructions for adjustment, maintenance and repair.


1. Fill the air bag orally until the cells are smooth.

2. After the cushion is inflated, the user can sit in a normal sitting position with correct arm and foot position.

3. Place a hand between the cushion and the user's tailbone to adjust the cushion into position.

4. Gently unscrew the valve by hand to deflate and adjust the pillow to the user.

#### WARNING:

In the event of a „serious incident” related to the product that has resulted, could have resulted, or is likely to result, directly or indirectly, in any of the following situations:

- a) the death of the patient, user or other person, or
- b) temporary or permanent deterioration in the health of the patient, user or other person, or
- c) a serious risk to public health

this „serious incident” must be reported to the manufacturer and to the competent authority of the Member State where the user or patient is present. For Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

#### WARNING:

In case of pain, allergic reactions or other disturbing and unclear symptoms related to the use of the drug, consult your doctor.

#### Indications

Hip ulcers and decubitus ulcers, rotation and other changes in the pelvis, hip arthroplasty, wheelchair patients.

#### Utilisation

Please dispose of municipal waste in accordance with the applicable disposal regulations.

#### CONTRAINDICATIONS

Physical or mental limitations (e.g. visual impairment) that prevent safe handling of the product.

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Model: .....  
Serial number: .....

Point of sale stamp

and a legible signature of the Seller: .....

1. ANTAR MEDICAL Sp. z o.o. ,03-068 Warsaw, ul. Zawiślńska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.
10. The warranty does not cover parts and components subject to natural wear and tear, such as wheels, tires, rubber pads, fuses, etc., the wear of which results from normal use of the product. The warranty period for batteries is 6 months.

ANTAR MEDICAL Sp. z o.o.  
03-068 Warsaw, ul. Zawiślńska 43  
Tel. 22 518 36 00, Fax 22 518 36 30  
www.antar.net

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