AT53032

INSTRUKCJA OBSŁUGI ORTEZA NADGARSTKA
USER MANUAL FOR WRIST ORTHOSIS
GEBRAUCHSANWEISUNG HANDGELENKSORTHESE
NÁVOD K POUŽITÍ ORTÉZY ZÁPĚSTÍ
NÁVOD NA POUŽITIE - ORTÉZA ZÁPÄSTIA
MODE D'EMPLOI DE L'ORTHÈSE DE POIGNET
GEBRUIKSAANWIJZING POLSORTHESE
INSTRUCCIONES DE USO ORTESIS DE MUÑECA
ISTRUZIONI PER L'USO DELL'ORTESI DI POLSO
BRUKSANVISNING FÖR HANDLEDSORTOS

IU_AT53032

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Thank you for purchasing our product, please read the instructions for use carefully. Do not use the product before reading and understanding this manual. If you do not understand the warnings, cautions and recommendations, please contact your health care professional or dealer to avoid damage to the product or bodily injury.

NOTE: Inspect all parts for damage caused in transit. If you notice such damage, DO NOT use the product. Contact the manufacturer for more information.

USE/INDICATIONS

The orthosis is indicated for all users who need help stabilizing the wrist.

Indicated for:

- carpal tunnel syndrome
- dislocation of the wrist joint or overuse syndrome
- after wrist surgery

CONTRAINDICATIONS

Physical or mental limitations (e.g., visual impairment) that prevent safe handling.

- The use of a wrist orthosis for open wounds or skin infections can worsen the skin condition, increase the risk of spreading infection, and impede healing.
- If the orthosis is too tight or worn for too long, it can restrict blood circulation to the wrist and hand area, which can lead to further complications such as pain, numbness or swelling.
- A brace that is not properly fitted can cause discomfort, skin abrasions or even clinical deterioration by improperly stabilizing the joint.
- People with sensory impairments in their hand may not feel pain or discomfort from improper use of the orthosis, which can lead to further tissue damage.
- In cases of advanced arthritis or other conditions that require greater mobility of the wrist for physiotherapy, an overly stiff orthosis may limit the necessary movements for rehabilitation.

FEATURES

- Cut wrist **orthosis** for left and right hand
- **Profiled** plastic insert and an elaborate fastening strap system provide complete support for the wrist and thumb
- made of breathable material warms the wrist and improves blood circulation

SIZE


The product comes in right and left hand versions

Measure the width of the metacarpophalangeal joint (cm): XS: 5.1–6.3 cm S:

6.3–7.4 cm M: 7.4–8.5 cm L: 8.5–9.5 cm XL: 9.5–10.7 cm

Wrist circumference (cm): XS: 10.2–12.4 cm S: 12.5–15.2 cm M: 15.2–20.3 cm

L: 17.8–20.3 cm XL: 20.3–22.8 cm

COMPONENTS:

Orthosis, Instructions for use

NOTES:

In the event of a product-related "serious incident" that directly or indirectly led to, may have led to, or may lead to any of the following:

- (a) death of a patient, user or other person or
- (b) temporary or permanent deterioration of the health of the patient, user or other person or
- (c) a serious threat to public health

the above "serious incident" should be reported to the manufacturer and to the competent authority of the Member State in which the user or patient resides. In the case of Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

NOTES:

If you experience pain, allergic reactions, or other disturbing symptoms that are unclear to you related to the use of the medical device, consult a healthcare professional.



NOTE: it is forbidden to use the product other than for its intended purpose .

NOTE: The manufacturer is not responsible for damage caused by negligent cleaning and maintenance, or resulting from failure to follow the recommendations in this manual.

TARGET PATIENT GROUP

Persons struggling with diseases, dysfunctions or injuries, for the treatment/rehabilitation or compensation of which this device is intended (see the section on the purpose of the device of this manual). The device may be purchased by the user independently or on the recommendation of a physician, therapist or other specialist. Whether purchasing the device on its own or as indicated by a doctor/therapist/other specialist, the available sizes/necessary functions and variants of the device, indications and contraindications for use, as well as information provided by the manufacturer, should be taken into account.

USE

<p>1. Adjust the appropriate size of the orthosis by measuring the width of the metacarpophalangeal joint. 2. Undo all the Velcro, insert your hand into the contoured rail as shown.</p>		<p>3. Fasten the Velcro securing the wrist as shown.</p>	
<p>4. Fasten the velcro stabilizing the fingers and make sure the hand is in the right position.</p>		<p>5. Fasten the thumb stabilizing velcro as shown. 6. Then fasten the Velcro stabilizing the forearm. 7. Check that there is adequate air circulation, if necessary, the Velcro fastening can be improved.</p>	

CLEANING AND CARE

Fasten all Velcro before washing. Hand wash in cold water using mild soap. Do not dry clean or tumble dry. Do not use bleach. Do not expose to direct sunlight.

STORAGE

Store the product in a dry cool room, away from moisture and direct sunlight. Not suitable for children under three years of age. Keep away from sources of fire.

METHOD OF DISPOSAL OF THE PRODUCT

Once the device is decommissioned, the medical device can be disposed of like regular municipal waste

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WARRANTY CARD

Model:
Serial number:

Point of sale stamp
and a legible signature of the Seller:

1. ANTAR MEDICAL Sp. z o.o., 03-068 Warsaw, ul. Zawiańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

 ANTAR MEDICAL Sp. z o.o.
03-068 Warsaw, ul. Zawiańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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