TENS KTR-2610

**INSTRUKCJA UŻYTKOWNIA, ELEKTROSTYMULATOR MIĘŚNI
USER MANUAL, MUSCLE ELECTROSTIMULATOR
BEDIENUNGSANLEITUNG, MUSKEL-ELEKTROSTIMULATOR
UŽIVATELSKÁ PŘÍRUČKA, SVALOVÝ ELEKTROSTIMULÁTOR
NÁVOD NA POUŽITIE, SVALOVÝ ELEKTROSTIMULÁTOR
NOTICE D'UTILISATION, ÉLECTROSTIMULATEUR MUSCULAIRE
GEBRUIKSAANWIJZING, SPIER-ELEKTROSTIMULATOR
MANUAL DE USUARIO, ELECTROESTIMULADOR MUSCULAR
MANUALE D'USO, ELETTROSTIMOLATORE MUSCOLARE
BRUKSANVISNING, MUSKELELEKTROSTIMULATOR**

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EN

Thank you for purchasing our product. Please read this user manual carefully. Do not use the device before you have read and fully understood this manual. If you do not understand the warnings, notices, or recommendations, please contact a healthcare professional or the seller to avoid damage to the product or personal injury. CAUTION: All parts must be checked for damage that may have occurred during transport. If any such damage is found, the product MUST NOT be used. For further information, please contact the seller.

TARGET PATIENT GROUP

People suffering from diseases, dysfunctions, or injuries for which this device is intended for treatment, rehabilitation, or compensation (see the section on intended use in this user manual). The product may be purchased independently by the user or based on the recommendation of a physician, therapist, or another specialist. Whether the product is purchased independently or recommended by a healthcare professional, available sizes, required functions and variants, indications and contraindications for use, as well as the information provided by the manufacturer must be taken into account.

PRODUCT INTRODUCTION

The muscle electrostimulator is a portable, battery-powered, multifunctional device.

The device offers TENS (Transcutaneous Electrical Nerve Stimulation), EMS (Electrical Muscle Stimulation), and FITNESS functions. It has 50 operating modes that can generate electrical impulses through electrodes placed on the skin, helping users benefit from body stimulation. The electronic module includes control elements such as an ON/OFF button, a display, a mode selection button, and intensity adjustment buttons.

The LCD screen displays the selected mode and program, intensity, stimulation frequency, and remaining operating time.

The device is equipped with accessories including electrodes, electrode wires, and batteries. The electrode wires are used to connect the electrodes to the main unit.

The electrodes comply with biocompatibility standards ISO 10993-5 (cytotoxicity) and ISO 10993-10 (irritation and sensitization) and are replaceable.

INTENDED USE OF THE PRODUCT

TENS (modes 1–20)

For temporary relief of pain associated with sore and aching muscles in the shoulders, waist, back, neck, arms, legs, and feet caused by physical exertion or daily household activities, by applying electrical current to stimulate nerves.


EMS (modes 21–40)

Intended for relaxation of muscle spasms, prevention or retardation of muscle atrophy, improvement of local blood circulation, and muscle re-education.

FITNESS (modes 41–50)

Improvement of abdominal muscle tone, strengthening of abdominal muscles, and shaping a firmer abdomen.

ELECTRODE PLACEMENT SUGGESTIONS

TENS MODE (1–20)

Recommended operating time for all TENS programs: 30 minutes.

Lower back pain 1 – electrode positions: 4–11

Lower back pain 2 – electrode positions: 4–11

Rheumatic pain 1 – electrode positions: 1–28, especially 27

Rheumatic pain 2 – electrode positions: 1–28, especially 27

Thigh pain 1 – electrode positions: 22–24

Thigh pain 2 – electrode positions: 22–24

Shoulder pain 1 – electrode positions: 1–9, 16–19

Foot pain 1 – electrode positions: 23–26, 28

Shoulder pain 2 – electrode positions: 1–9, 16–19

Rheumatic pain 3 – electrode positions: 1–28, especially 27

Shoulder pain 3 – electrode positions: 1–9, 16–19

Thigh pain 3 – electrode positions: 22–24

Foot pain 2 – electrode positions: 23–26, 28

Shoulder pain 4 – electrode positions: 1–9, 16–19

Lower back pain 3 – electrode positions: 4–11

Foot pain 3 – electrode positions: 23–26, 28

Lower back pain 4 – electrode positions: 4–11

Thigh pain 4 – electrode positions: 22–24

Thigh pain 5 – electrode positions: 22–24

Shoulder pain 5 – electrode positions: 1–9, 16–19

EMS mode (21–40)

Application: muscle relaxation and support of physical recovery.

Possible electrode positions: 1–28.

Recommended operating time: 30 minutes.

FITNESS mode (41–50)

Application: strengthening of local muscle groups and development of muscle strength.

Possible electrode positions: 1–28, especially 12–14 and 16–24.

Recommended operating time: 15 minutes.

WARNING

- a) The long-term effects of chronic electrical stimulation on the human body are unknown.
- b) Patients with a cardiac pacemaker – use is prohibited.
- c) Patients with suspected or confirmed epilepsy, heart disease, pregnant women, or women during menstruation – use is prohibited.
- d) Do not use on areas with bleeding, recent injuries or fractures, during wound healing, or after surgery.
- e) Do not use on patients insensitive to heat or electrical skin stimulation.
- f) Do not use on persons with impaired consciousness or on children with disabilities.
- g) Do not use in cases of swelling, skin infections, skin diseases, or phlebitis.
- h) Do not use on wet skin after bathing, during heavy sweating, or before sleep.
- i) Do not use while driving or during movement.
- j) Do not use on persons with metal allergies.


ADVERSE REACTIONS

Skin irritation or a burning sensation may occur under the electrodes during stimulation.

Headaches or other painful sensations may occur during or after electrical stimulation near the eyes, head, or face. If adverse reactions occur, stop using the device immediately and consult a physician.

MAIN PRODUCT STRUCTURE

1. Functional modules
2. Channel 1 intensity
3. Program selection (up)
4. Channel 1 intensity adjustment
5. Functional module switch / stimulation time setting
6. Program
7. Electrical pulse level
8. Frequency
9. Channel 2 intensity
10. Operating time
11. Program selection (down)
12. Channel 2 intensity adjustment
13. On/Off button

MAIN PRODUCT FEATURES

The product features an ergonomic design combined with modern, aesthetic styling and a compact size for easy portability. Low-frequency electrical impulses generate an effective composite energy field.

The device offers multifunctional operation and allows simultaneous stimulation of multiple areas.

Multiple functional modules are available, with adjustable therapy duration and adjustable intensity.

Equipped with an LCD display.

LIST OF COMPONENTS: Main unit, charging cable, electrodes – 2 pairs, electrode wires – 2 pcs, user manual.

PRODUCT REQUIREMENTS AND MAIN PARAMETERS

Power Supply

Power supply: lithium battery 3.7 V / 250 mAh

Device safety classification: BF class, BF type

Main Technical Parameters

a) Impulse frequency:

TENS: 20–100 Hz

EMS: 1–15 Hz

FITNESS: 2–16 Hz

b) Pulse width:

TENS: 120 μ s

EMS: 200 μ s

FITNESS: 200 μ s

c) Pulse waveform: square wave

d) Impulse charge at maximum output amplitude: $\geq 7 \mu$ C

e) Maximum energy per single pulse: ≤ 300 mJ

f) Influence of open circuit and short circuit on output:

The device is designed to operate under both open-circuit and short-circuit conditions without degradation of performance.

g) Output amplitude adjustment: 0–16 levels

h) Therapy timer: 15 minutes, 30 minutes, 45 minutes ($\pm 10\%$ tolerance)

i) Safety classification: internal electric power source, BF-type equipment

j) Main unit dimensions: 122 mm \times 55 mm \times 15.7 mm

Electrodes:

Size: 50 \times 50 mm

Effective area per pad: 25 cm²

Material: hydrogel

k) Service life:

Device shelf life: 2 years

Small electrode pad lifespan: 80 uses

Electrode shelf life: 2 years

ENVIRONMENTAL REQUIREMENTS

(Transport and storage conditions between uses)

a) Normal operating conditions:

Ambient temperature: +5°C to +40°C

Relative humidity: 15%–93% RH

Atmospheric pressure: 700 hPa–1060 hPa

b) Storage conditions:

Temperature: –25°C to +55°C

Humidity: 0%–93% RH

Atmospheric pressure: 700 hPa–1060 hPa

c) Transport conditions:

Temperature: –25°C to +55°C

Humidity: 0%–93% RH

Atmospheric pressure: 700 hPa–1060 hPa

Note:

If the product has been stored at minimum or maximum temperatures, allow it to remain at normal ambient temperature (approx. 20°C) for about 2 hours before use to restore proper functionality.

POWER OPERATION

Charge the lithium battery by connecting the device to the charger.

Connect all components of the set.

Connect the main unit to the electrodes using the electrode wires.

Power On / Off

Press the On/Off button to turn on the device. The power indicator will light up, the operating program will start, and the LCD backlight will activate.

Press the On/Off button again to turn off the device and switch off the display.

FUNCTIONAL OPERATION

Attach the electrodes to the areas intended for stimulation. If the skin is sweaty, clean it before use.

Press the On/Off button to turn on the main unit. The device defaults to TENS mode (20 preset programs).

Use the ▲ and ▼ buttons to select a program:

▲ – program up

▼ – program down

Use the + and – buttons to adjust intensity for channels 1 and 2.

Each channel has 16 intensity levels:

+ increases intensity

– decreases intensity

Mode Switching

Short press the F button to switch between TENS, EMS, and FITNESS modes.

Therapy Time Setting

Long press the F button (approx. 1.5 seconds) to set the stimulation time: 15 minutes, 30 minutes, or 45 minutes.

BATTERY USE

The device is powered by a rechargeable lithium battery.

If the LCD brightness decreases significantly or intensity changes become weak during button operation, the battery level is low and must be recharged immediately.

If the device has not been used for more than one month, fully charge it before reuse.

Do not use the device in environments above 45°C, as this may affect performance and shorten battery life. Dispose of used batteries in accordance with local environmental regulations.

SAFETY PRECAUTIONS

The therapeutic device must not be used simultaneously with high-frequency (HF) equipment to avoid burns or equipment damage.

Simultaneous use with HF equipment may cause burns at the electrode or massage area and damage the device.

Operation within 1 meter of shortwave or microwave equipment may result in unstable performance.

Applying electrodes near the chest may increase the risk of cardiac fibrillation.

Do not modify the device without manufacturer authorization.

Battery replacement must be performed only by authorized service personnel.

The device is not intended for use by persons unable to care for themselves, mentally impaired persons, or unconscious individuals.

Simultaneous connection to high-frequency surgical equipment may cause burns and device damage.

Do not apply electrodes on or near the chest, eyes, head, or front of the neck (especially the carotid sinus).

Prevent inhalation or accidental swallowing of small parts (including AAA/LR03 batteries).

Do not use accessories not approved by the manufacturer.

CONTRAINDICATIONS

Pregnant women, women during menstruation, persons with sensitive skin, heart disease, abnormal blood pressure, malignant tumors, cerebrovascular disease, acute illness, or those under medical treatment should consult a physician before use.

The product is contraindicated for persons with impaired skin sensation or insensitivity to heat.

Do not use during bathing, heavy sweating, or sleep.

In patients with cerebral hemorrhage, use is prohibited during unstable phases; post-stroke patients may use only under medical supervision.

The product is contraindicated in cases of purulent inflammation, acute blood poisoning, persistent high fever, and acute cardiovascular or cerebrovascular diseases.

If discomfort or skin reactions occur, discontinue use immediately and consult a physician.

PRECAUTIONS

Do not apply the device near the heart, head, eyes, front of the neck (especially carotid artery), lower back, oral cavity, or diseased skin areas.

Switch the device off and restart it after changing therapeutic settings or repositioning to prevent excessive stimulation.

Do not use on children or persons unable to communicate sensations clearly.

Disconnect the power plug when not in use.

Do not use simultaneously with other medical electronic devices such as pacemakers, artificial heart-lung systems, ECG equipment, or similar devices, as this may be life-threatening.

ADDITIONAL WARNINGS

Do not use the product in locations exposed to direct sunlight, high temperatures, flammable materials, electromagnetic radiation, or moisture.

Do not disassemble, repair, or modify the therapeutic device, as this may cause malfunction or electric shock.

Place the device in a location that allows easy access to the power plug so that the power supply can be quickly disconnected in an emergency.

Before each use, inspect the device and avoid using damaged cables or components.

If the device performance decreases, clean it with a dry, soft cloth.

Before each use, check the electrodes—dirty electrodes may reduce effectiveness or cause operational issues.

The product contains connecting wires; prevent improper use by children to avoid the risk of strangulation.

NOTICE:

If a “serious incident” related to the device occurs that directly or indirectly has led, could have led, or may lead to any of the following events:

- a) death of the patient, user, or another person, or
- b) temporary or permanent deterioration of the health of the patient, user, or another person, or
- c) a serious threat to public health,

the incident must be reported to the manufacturer and to the competent authority of the Member State in which the user or patient resides.

NOTICE: If pain, allergic reactions, or other alarming or unclear symptoms related to the use of the medical device occur, consult a healthcare professional.

NOTICE: Use of the product for purposes other than its intended use is prohibited.

PRODUCT MAINTENANCE

Before switching on the device, ensure that the battery is charged.

If the device does not start properly, switch it off and turn it on again after a few seconds. If the problem persists, check the fuse.

If the device operates but no stimulation is felt, check whether the electrodes are in proper contact with the skin and whether there are any obstructions (hair, clothing, etc.).

Before cleaning, switch off the device and disconnect the electrodes from the power supply. Clean with a soft cloth slightly moistened with water, then dry thoroughly.

After each use, clean the electrodes. In case of heavy contamination, a soft cloth with a small amount of medical alcohol (75% concentration) may be used.

Store the device in a dry, well-ventilated place away from moisture and heat sources.

Handle the device carefully during transport and avoid shocks or impacts.

Regularly check the battery condition and overall device performance. If abnormalities occur, exercise caution and promptly submit the device for servicing.

The device housing should be adequately protected against mechanical damage, especially abrasion and scratches.

Persons without appropriate qualifications must not disassemble the main unit. Doing so may result in electric shock, device damage, or an accident posing a risk to health or life.

PRODUCT DISPOSAL

The therapeutic device is a medical device. If it becomes worn out or is no longer suitable for use, it must be disposed of in accordance with applicable local regulations concerning the disposal of medical devices and electrical waste.

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Shenzhen Kentro Medical Electronics Co.,Ltd 2nd Floor, No. 11, Shanzhuang Road, Xikeng Village, Yuanshan Street, Longgang District, Shenzhen City, Guangdong Province, China
[Tel]: +86(755)33825998, [Fax]: +86(755)33825996

Share Info GmbH

Am Schulzentrum 12, 41564 Kaarst, Germany
Tel: 0179 5666 508 Email: eu-rep@share-info.com

IMPORTER: ANTAR MEDICAL Sp. z o.o.
03-068 Warszawa, ul. Zawisłańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net
e-mail: antar@antar.net

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v1-02.02.2026

WARRANTY CARD

Model:
Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. ANTAR MEDICAL Sp. z o.o. ,03-068 Warsaw, ul. Zawiślńska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.
10. The warranty does not cover parts and components subject to natural wear and tear, such as wheels, tires, rubber pads, fuses, etc., the wear of which results from normal use of the product. The warranty period for batteries is 6 months.

 Shenzhen Kentro Medical Electronics Co.,Ltd 2nd Floor, No. 11, Shanzhuang Road, Xikeng Village, Yuanshan Street, Longgang District, Shenzhen City, Guangdong Province, China
[Tel]: +86(755)33825998, [Fax]: +86(755)33825996
Share Info GmbH
Am Schulzentrum 12, 41564 Kaarst, Germany
Tel: 0179 5666 508 Email: eu-rep@share-info.com

 ANTAR MEDICAL Sp. z o.o.
03-068 Warsaw, ul. Zawiślńska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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