**AT52337**

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En

Thank you for purchasing our product. Please read this instruction manual carefully before use. Do not use the product before reading and fully understanding this manual. If you do not understand the warnings, notes, and recommendations, contact a healthcare professional or the seller to avoid product damage or personal injury. **WARNING:** All parts must be checked for damage that may have occurred during transport. If any such damage is found, **DO NOT** use the product. Contact the seller for further information.

#### **INTENDED PATIENT GROUP**

Individuals suffering from diseases, dysfunctions, or injuries for the treatment, rehabilitation, or compensation of which this product is intended (see the section on the intended use of the product in this user manual).

The product may be purchased independently by the user or based on the recommendation of a physician, therapist, or other specialist. Whether the product is purchased independently or recommended by a physician/therapist/other specialist, available sizes, required functions and variants of the product, indications and contraindications for use, as well as the information provided by the manufacturer must be taken into account.

**INTENDED USE:** An electric wheelchair is a vehicle structurally designed for the mobility of a person with a disability. It is powered by a motor whose design limits the maximum speed to walking pace.

#### **INDICATIONS**

**Severe mobility limitations:** Individuals with significant limitations in muscle strength, for example due to neurological or muscular diseases or pain syndromes, which prevent or significantly hinder independent mobility.

**Neurodegenerative diseases:** Individuals suffering from conditions such as multiple sclerosis, Parkinson's disease, muscular dystrophies, or ALS, which lead to progressive loss of muscle control and strength.

**Fatigue:** Individuals for whom prolonged walking or manual operation of a wheelchair is excessively demanding due to illness or health condition, causing rapid fatigue and exhaustion.

**Age-related limitations:** Older people who, due to age-related reductions in mobility and strength, may benefit from an electric wheelchair to maintain independence and mobility.

**Rehabilitation and recovery periods:** Individuals undergoing rehabilitation after injuries or surgeries who are temporarily unable to move independently may use an electric wheelchair as part of the treatment and recovery process.

#### **CONTRAINDICATIONS**

Physical or mental limitations (e.g. visual impairment) that prevent safe handling of the product.

Individuals who have difficulty judging distances or spatial relationships may have problems safely maneuvering an electric wheelchair.

Electric wheelchairs require a certain level of manual ability and coordination to operate the joystick or other control devices. Persons who are unable to operate these controls effectively may not be able to use the electric wheelchair safely.

Although many individuals with cognitive impairments can safely use electric wheelchairs, in some cases severe impairments may prevent safe operation. Individuals who are unable to maintain a stable seated position may require additional support systems to safely use an electric wheelchair. Epilepsy or other conditions that may lead to sudden loss of consciousness or reaction may constitute a contraindication. If a person has impaired motor coordination or is unable to independently operate the wheelchair controls, this may be a contraindication. Individuals with disorders of consciousness that may affect their ability to safely operate an electric wheelchair may also be unsuitable for its use. In some cases, serious heart conditions or blood pressure disorders may limit the safe use of an electric wheelchair, especially if emotions or stress may trigger adverse reactions. Conditions such as advanced multiple sclerosis or other diseases affecting muscle control may require specialized adaptive solutions in wheelchairs. Problems with depth perception, spatial perception, or general orientation may significantly reduce the safety of using an electric wheelchair. Generally poor health or recent surgeries, especially involving the spine or lower limbs, may require postponing the use of the wheelchair until health has improved.

**LIST OF COMPONENTS:** Wheelchair, battery, charger, user manual

1. Push handle
2. Armrest
3. Motor
4. Drive wheel
5. Controller
6. Battery
7. Footrest

#### **TECHNICAL SPECIFICATIONS**

##### **CONSTRUCTION**

Frame: aluminum alloy

Seat, backrest and cushion material: foam, polyester mesh fabric

##### **OVERALL DIMENSIONS**

Length: 104 cm

Width: 65 cm

Height: 98 cm

##### **FOLDED DIMENSIONS**

Folded length: 70 cm

Folded width: 45 cm

Folded height: 75 cm

#### SEAT

Seat width: 46 cm

Seat depth: 45 cm

Seat thickness: 7 cm

Seat height from the ground (front): 45 cm

Seat tilt angle: 3°

#### BACKREST AND ARMRESTS

Backrest height: 57 cm

Backrest tilt angle: 13.5°

Armrest height from seat: 31 cm

Armrest height from ground: 76 cm

#### WHEELS

Front wheels: solid, 8", PU

Rear wheels: pneumatic, 12", rubber

Tire pressure: 2.0

#### DRIVE AND ELECTRONICS

Motor: 2 × 250 W, 24 V, brushed

Battery: lithium, 24 V / 12 Ah

Battery weight: 2.5 kg

Charger: input AC 220 V, 50 Hz / output 24 V, 2 A

Maximum controller output current: 1.5 A

#### PERFORMANCE PARAMETERS

Maximum speed: 6 km/h

Maximum range per charge: 20 km

Braking distance at 6 km/h: ≤ 1200 mm

Minimum turning radius: 900 mm

Maximum climbing ability: ≤ 13°

Maximum braking ability on an incline: ≤ 13°

#### LOAD AND WEIGHT

Maximum user weight: 120 kg

Wheelchair weight with battery: 29 kg

#### ADDITIONAL INFORMATION

Can the wheelchair be used as a seat in a motor vehicle: No

NOTE: Individual product dimensions may vary by ±2 cm.

This symbol indicates the maximum user weight.

#### ASSEMBLY

Remove the electric wheelchair and place it on a flat surface. With one hand, hold the backrest tube, and with the other hand press down on the seat, unfolding the wheelchair firmly until it is fully open. Make sure that the locking latch is properly engaged.

Controller installation: Insert one end of the controller extension tube into the opening in the tube under the wheelchair armrest. Then tighten the knob and align the power cable in a straight line.

Adjustment of anti-tip Wheel: Slide the anti-tip wheel tube backward to more effectively prevent the wheelchair from tipping backward.

Backrest height adjustment: Loosen the two knobs located at the rear of the backrest and adjust the backrest height to the desired position.

Clutch operation

In the event of an electric wheelchair malfunction or battery discharge, switch the wheelchair to manual mode and move it to a safe location.

Switching from electric to manual mode: move the clutch lever to the manual position.

Switching from manual to electric mode: move the clutch lever to the electric position.

Opening the armrest: Press and hold the plastic button on the side of the armrest, grasp the armrest tube, and lift it upward to make getting in and out of the wheelchair easier. Before driving, make sure the armrest is locked.

Use of safety Belt: Always fasten the safety belt when using the electric wheelchair.

#### **CHARGER AND BATTERY**

The charger is used to charge the internal power source of the wheelchair. The wheelchair must not be used while charging.

Charger parameters: Input voltage: AC 220 V, 50 Hz/ Output voltage: DC 24 V/ Output current: 2 A/ Ingress protection rating: IPX3

Charger use:

Connect the charger plug to the socket located under the controller.

Charging procedure: Check that the charger socket is not blocked. Make sure the controller is switched off. Connect the charger output plug to the socket under the controller, then connect the input plug to a 220 V household power outlet – charging will begin.

When all green indicator lights on the charger are illuminated, charging is complete.

After charging, disconnect both the input and output plugs of the charger.

During charging, the charger should first be connected to the battery and only then to the power outlet. The maximum charging time is 6–8 hours. When the green indicator light is on, charging is complete. After disconnecting the charger from the power supply, also disconnect the plug on the battery side to prevent battery discharge. If the wheelchair is not used for a long period, it is recommended to charge the battery once a month to extend its service life.

## BATTERY REMOVAL

Disconnect the cable connector from the battery box. Press the button located on the upper front part of the battery box and slide it out. To install the battery, perform the above steps in reverse order.

## NOTE

After charging is completed, the charging indicator turns green. Do not interrupt the charging process before it is finished.

After charging is completed (usually after 6–8 hours), switch off the charger power supply and unplug it. If the cable is not disconnected, the battery will slowly discharge. The maximum charging time must not exceed 12 hours, as there is a risk of overcharging.

If the wheelchair is not used for a long period of time, the battery should be charged once a month to keep it fully charged.

Charging time depends on the ambient temperature; in winter, charging takes longer.

Note – safety rules during charging

To avoid hazards during charging, observe the following rules:

Use only the standard charger supplied by the manufacturer. The use of non-standard chargers may be hazardous. Modifying or altering charging cables is prohibited. The manufacturer is not responsible for accidents or damage resulting from such actions.

The charging area should be well ventilated and protected from direct sunlight and moisture.

Do not cover the charger or the wheelchair with fabrics or other objects during charging.

Do not charge the battery at temperatures below –10°C or above 50°C, as the charger may not operate correctly and the battery may be damaged.

The charger heats up during charging – this is a normal phenomenon related to heat dissipation.

The charger is not waterproof.

During charging, do not place the charger on flammable materials such as oils, pedals, or seats.

## WARNING

During charging, keep the wheelchair away from high temperatures and open flames, as this may cause ignition or explosion of the battery. Do not charge the battery in enclosed spaces.

During charging, the wheelchair should be placed in an open area, away from flammable materials, to ensure fire safety.

## SAFE USE AND MAINTENANCE OF THE BATTERY

Incorrect battery replacement poses a risk of explosion. When replacing the battery, use only a battery of the same type as the original or one recommended by the manufacturer, and ensure correct connection of the positive and negative terminals.

To extend battery life, charge it fully on a regular basis. If the wheelchair is not used for a short period, the battery should be fully charged and ready for use. In the case of long-term non-use, it is recommended to charge the battery once a month.

## CONTROLLER

1. Joystick
2. Battery level indicator
3. Power on/off button
4. Speed indicator
5. Acceleration button
6. Mute button
7. Horn
8. SOS button

After switching on the power, the power indicator lights up, the battery indicator lights up sequentially from left to right, and the buzzer emits a sound signal – the controller is operating correctly. After switching off the power, all LEDs turn off and the controller does not generate an output signal.

Warning – battery indicator

The battery indicator has three colors: red, yellow, and green. All five LEDs lit indicate a fully charged battery. If only the red LED is lit, this indicates a low battery level – stop using the wheelchair immediately and charge the battery.

If the error persists after releasing the joystick, this may indicate damage to the joystick. In such a case, do not use the wheelchair and contact service.

Sleep mode: If the joystick is not used for 10 minutes, the controller automatically switches off the power and enters sleep mode. Pressing the power button wakes up the control system. In sleep mode, power consumption is minimal and does not damage the battery.

Speed adjustment: The driving speed can be adjusted according to user preference and environmental conditions. The controller offers 5 speed levels – level 5 is the fastest and level 1 the slowest. Speed can be adjusted using the acceleration and deceleration buttons.

Horn operations: Pressing the horn button produces a continuous sound signal, which stops when the button is released.

The battery level indicator shows that the wheelchair power is on and provides approximate information about the remaining energy level.

All LEDs lit – battery fully charged

Yellow or red LEDs lit – battery requires charging

Only red LED lit – very low energy level, charge the battery immediately

Joystick operation

Joystick movement determines the direction of wheelchair travel:

Forward – move forward/Backward – move backward/ Left – rotate left on the spot/ Right – rotate right on the spot/The degree of joystick deflection determines the driving speed in the given direction.

#### IMPORTANT NOTICE – Serious incidents

In the event of a “serious incident” related to the medical device that directly or indirectly led, could have led, or may lead to any of the following events:

- a) death of a patient, user, or other person, or
- b) temporary or permanent deterioration of the health of a patient, user, or other person, or
- c) a serious public health threat,

the above “serious incident” must be reported to the manufacturer and to the competent authority of the Member State in which the user or patient resides.

NOTE: If pain, allergic reactions, or other disturbing or unclear symptoms related to the use of the medical device occur, consult a healthcare professional.

NOTE: It is prohibited to use the product in any manner other than its intended purpose.

#### CLEANING AND MAINTENANCE

User-performed maintenance includes cleaning, checking the tires, the battery, and the charging process. More advanced maintenance should be carried out by service personnel. It is recommended to inspect the wheelchair every two months and perform a full technical inspection once a year.

Cleaning the wheelchair

Regularly remove dust from the wheelchair surfaces. Use wax polish or a clean, soft cloth for cleaning.

Tires

Regularly check the wear of pneumatic tires and replace them when the tread depth reaches 1 mm.

Battery

The battery should be charged regularly. Fully discharging the battery before recharging is not recommended, as it shortens battery life.

#### STERILIZATION AND DISINFECTION

The wheelchair is not a sterile medical device. During daily use, it should be kept clean.

If the wheelchair becomes contaminated with blood, secretions, excretions, or other substances, wipe it with a disinfectant based on ammonium salts at a concentration of 1000–2000 mg/L, leave it for 30 minutes, then wipe with a clean, dry cloth and allow it to dry.

#### WHEELCHAIR DAMAGE AND SELF-CHECK

If any abnormalities occur during use, switch off the power and check the wheelchair.

Symptom: complete loss of power and the controller indicators are off.

Step 1: Check whether the controller plug is loose or disconnected.

Step 2: Check the connection between the battery box and the controller. If necessary, reconnect the plug by holding the connector, not the cable, to avoid damage to the wiring.

If power is not restored after performing the above steps or if other concerns arise, contact service immediately.

#### STORAGE AND TRANSPORT CONDITIONS

The electric wheelchair should be stored and transported under the following conditions:

Temperature range: 5°C – 40°C

Relative humidity: ≤ 80%

Atmospheric pressure: 860 hPa – 1060 hPa

Rated voltage

Charger operating voltage: AC 220 V, 50 Hz

Internal power supply voltage: DC 24 V

#### DISPOSAL

This product contains electronic components. At the end of its service life, it must not be disposed of with household waste. For proper disposal, return the product to an appropriate collection point for electrical and electronic equipment. Proper disposal helps prevent potential harm to the environment and human health in accordance with applicable legal regulations on the recycling of electrical and electronic equipment.

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Model: .....  
Serial number: .....

Point of sale stamp

and a legible signature of the Seller: .....

1. ANTAR MEDICAL Sp z o.o, 03-068 Warsaw, ul. Zawisłańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.
10. The warranty does not cover parts and components subject to natural wear and tear, such as wheels, tires, rubber pads, fuses, etc., the wear of which results from normal use of the product. The warranty period for batteries is 6 months.

ANTAR MEDICAL Sp. z o.o  
03-068 Warsaw, ul. Zawisłańska 43  
Tel. 22 518 36 00, Fax 22 518 36 30  
www.antar.net

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