AT51133

INSTRUKCJA UŻYTKOWNIA, KULA ŁOKCIOWA
USER MANUAL, ELBOW CRUTCH
GEBRAUCHSANWEISUNG, UNTERARMGEHSTÜTZE
NÁVOD K POUŽITÍ, FRANCOUZSKÁ HŮL
NÁVOD NA POUŽITIE, FRANCÚZSKA BARLA
MODE D'EMPLOI, BÉQUILLE D'AVANT-BRAS
GEBRUIKSAANWIJZING, ELLEBOOGKRUK
MANUAL DE USO, MULETA DE ANTEBRAZO
ISTRUZIONI PER L'USO, STAMPELLA ANTIBRACHIALE
BRUKSANVISNING, UNDERARMSKRYCKA

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En Thank you for purchasing our product. Please carefully read the user manual before use. Do not use the product before reading and understanding this manual. If you do not understand the warnings, precautions, and recommendations, please consult a healthcare professional or the seller to avoid product damage or personal injury.

WARNING: Check all parts for damage that may have occurred during transport. If such damage is found, DO NOT use the product. For more information, contact the seller.

TARGET GROUP OF PATIENTS

Individuals with diseases, dysfunctions, or injuries for which this product is intended as part of treatment, rehabilitation, or compensation (see the section on intended use in this manual). The product may be purchased independently by the user or based on recommendations from a doctor, therapist, or another specialist. In both cases, it is important to consider the available sizes/necessary functions and product variations, the indications and contraindications for use, as well as the information provided by the manufacturer.

INTENDED USE: The elbow crutch is designed to assist individuals with reduced mobility in walking. It facilitates movement, improves stability, and reduces the load on the lower limbs and spine. It can be used both at home and in healthcare facilities.

INDICATIONS: Individuals who have difficulty walking independently. Patients undergoing rehabilitation after lower limb injuries Or orthopedic surgery. Elderly individuals requiring additional support for balance. Patients with chronic conditions that limit mobility.

CONTRAINDICATIONS

Physical or mental limitations (e.g., vision impairment) preventing safe use of the product. Patients unable to safely support themselves on the crutch (e.g., very weak muscles, lack of arm strength). Individuals with balance disorders requiring specialized mobility aids. Do not use if the crutch is damaged (cracks, loose parts, worn rubber tip). The product is not intended for climbing, overcoming obstacles, or as a support for sitting.

LIST OF COMPONENTS: Elbow crutch, User manual

PRODUCT FEATURES

Lightweight and durable aluminum construction
Adjustable height (10 levels of adjustment, range 92.5 – 115 cm)
Ergonomic plastic handle for user comfort
Non-slip rubber tip for stability
Foldable design for easy transport and storage

TECHNICAL SPECIFICATIONS

Height adjustment range: 92.5 – 115 cm
Total width: 12 cm
Total length: 17.5 cm
Net weight: 0.48 kg
Tube diameter: upper – 22 mm, lower – 19 mm
Maximum user weight: 100 kg

This symbol indicates the maximum user weight.

ASSEMBLY METHOD

Open the package and remove the parts. Press the locking ball and attach the joint of the crutch.

Press the locking ball and loosen the rotating clamp; adjust the crutch to the desired height.

Release the locking ball, tighten the rotating clamp, and secure the height.

To fold, press the locking ball at the joint and collapse the crutch.

INSTRUCTIONS FOR USE

Adjust the crutch to the appropriate height – the arm should be slightly bent at the elbow (approx. 15–20°) when holding the handle. Lock the chosen height with the adjustment mechanism (by pressing the adjustment buttons). Place the forearm in the cuff and grip the handle with your hand. While walking, move the crutch alternately with the opposite leg (e.g., right crutch – left leg). Ensure the rubber tip makes full contact with the ground and is not worn out. After use, the crutch can be folded for easier storage and transport.

NOTES

Before use, ensure all parts are in good condition, connections are stable, and adjustment mechanisms are properly secured.

Make sure the locking ball fully snaps into the adjustment hole – only then is the height stable.

Use with caution on slippery or inclined surfaces to avoid slipping.

The product is intended solely to support walking – do not use it for other purposes.

Do not use as a toy or as a support for weights.

WARNING

If a “serious incident” related to the product occurs that directly or indirectly has led, could have led, or may lead to any of the following:

- a) death of the patient, user, or another person, or
- b) temporary or permanent serious deterioration of health of the patient, user, or another person, or
- c) a serious public health risk,

this “serious incident” must be reported to the manufacturer and the competent authority of the Member State where the user or patient resides. If pain, allergic reactions, or other concerning or unclear symptoms occur while using the product, consult a healthcare professional. Do not use the product in any way other than its intended purpose.

CLEANING: Clean with a soft, damp cloth and mild detergent. Avoid using corrosive agents, bleach, or solvents. Regularly check the condition of the rubber tip – replace it if worn. Check the tightness of the fittings and tighten screws if necessary.

STORAGE: Store in a dry place and avoid direct sunlight.

DISPOSAL: Once the product is no longer in use, dispose of it as regular household waste. Follow local waste segregation regulations. Protect the environment.

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Model:
Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. ANTAR MEDICAL Sp. z o.o., 03-068 Warsaw, ul. Zawisłańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR MEDICAL Sp. z o.o.
03-068 Warsaw, ul. Zawisłańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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