AT04003(K)

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EN

Thank you for purchasing our product. Please read these instructions for use carefully. Do not use the product before reading and understanding these instructions. If you do not understand the warnings, cautions, and recommendations, contact a healthcare professional or the seller to avoid product damage or personal injury. CAUTION: Check all parts for damage. Before first use, carefully inspect the technical condition of the product and ensure that all parts are complete and properly assembled (if applicable, pay particular attention to screws, nuts, rivets, and other fastening and structural elements). Make sure the product is not defective or damaged. If any damage or missing parts are noticed, DO NOT use the product. Contact the seller for further information.

INTENDED PATIENT GROUP

Persons suffering from diseases, dysfunctions, or injuries for whose treatment, rehabilitation, or compensation this product is intended (see the intended purpose section of these instructions for use). The product may be purchased independently by the user or based on the recommendation of a doctor, therapist, or another specialist. Whether purchased independently or recommended by a healthcare professional, the available sizes, required functions and variants of the product, indications and contraindications for use, as well as the information provided by the manufacturer should be taken into consideration.

INTENDED PURPOSE: The product is intended for stabilization and support relief of the upper limb.

INDICATIONS

Indications for use include, among others:

mild sprains,
elbow joint dislocations,
post-traumatic rehabilitation period,
support of plaster casts.

The product allows the limb to remain in a resting position while limiting its mobility according to medical recommendations.

The product is intended for children over 3 years of age.

CONTRAINDICATIONS

Do not use the product:

in the case of open wounds or skin lesions in the area of contact with the product,
in the case of allergy or hypersensitivity to the materials used in the product,
without consulting a doctor in the case of serious injuries requiring specialist treatment,
in a manner inconsistent with its intended purpose or medical staff recommendations,
in children under 3 years of age.

CAUTION: The product contains small parts. Choking hazard. The product should only be used under adult supervision.

LIST OF COMPONENTS: Sling, instruction manual

PRODUCT FEATURES

Base material: spacer fabric
Suitable for left and right arm
Velcro detachable panda head
Movable eyes attached with Velcro
4 sizes: XXXS, XXS, XS, S

SIZES

To select the correct size, measure the length of the forearm:

XXXS: (17-19cm) XXS: (19 - 20.7 cm) XS: (21 - 23 cm) S: (23.5 - 25.4 cm)

MATERIAL COMPOSITION: POLYESTER – 73% POLYAMIDE – 15% POLYURETHANE – 7% COTTON – 3% POLYSTYRENE – 2%

USE

Insert the injured limb into the sling and place the supporting strap over the head, positioning it on the healthy shoulder. Adjust the strap length so that the arm is bent at the elbow joint.

CAUTION:

In the event of a “serious incident” related to the product that directly or indirectly led, could have led, or may lead to any of the following:

- a) death of the patient, user, or another person, or
- b) temporary or permanent deterioration of the health of the patient, user, or another person, or
- c) a serious public health threat,

the above-mentioned “serious incident” should be reported to the manufacturer and the competent authority of the Member State in which the user or patient resides. In the case of Germany, the competent authority is the Federal Institute for Drugs and Medical Devices (BfArM). CAUTION: If pain, allergic reactions, or other disturbing symptoms related to the use of the medical device occur, consult a healthcare professional. CAUTION: It is prohibited to use the product in any way other than its intended purpose.

CLEANING AND MAINTENANCE

The product should be hand washed. Recommended hand washing temperature: 30–40°C. Do not bleach, do not use strong detergents or chemical agents. Do not tumble dry, iron, or dry clean. Dry flat in a spread-out position and gently squeeze out excess water before drying.

STORAGE: Store the product in a dry place. Do not expose it to direct sunlight.

DISPOSAL: After the product is withdrawn from use, dispose of it as regular household waste. Follow local waste segregation regulations.

Protect the environment.

SYMBOL EXPLANATION:

Product is flammable – keep away from open flames.

Not suitable for children aged 0–3 years.

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Model:
Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. ANTAR MEDICAL Sp. z o.o. ,03-068 Warsaw, ul. Zawisłańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.
10. The warranty does not cover parts and components subject to natural wear and tear, such as wheels, tires, rubber pads, fuses, etc., the wear of which results from normal use of the product. The warranty period for batteries is 6 months.

ANTAR MEDICAL Sp. z o.o.
03-068 Warsaw, ul. Zawisłańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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