AT53116_AT04303

INSTRUKCJA UŻYTKOWANIA, ORTEZA NADGARSTKA I KCIUKA
INSTRUCTIONS FOR USE, WRIST AND THUMB ORTHOSIS
GEBRAUCHSANWEISUNG, HANDGELENK- UND DAUMENORTHESE
NÁVOD K POUŽITÍ, ORTÉZA ZÁPĚSTÍ A PALCE
NÁVOD NA POUŽITIE, ORTÉZA ZÁPÄSTIA A PALCA
MODE D'EMPLOI, ORTHÈSE DU POIGNET ET DU POUCE
GEBRUIKSAANWIJZING, POLS- EN DUIMORTHESE
INSTRUCCIONES DE USO, ORTESIS DE MUÑECA Y PULGAR
ISTRUZIONI PER L'USO, ORTESI DEL POLSO E DEL POLLICE
BRUKSANVISNING, HANDLED- OCH TUMORTOS

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EN

Thank you for purchasing our product. Please read the instructions for use carefully. Do not use the device before reading and understanding these instructions. If you do not understand the warnings, precautions and recommendations, contact a healthcare professional or the seller to avoid damage to the product or personal injury. **CAUTION:** Check all parts for damage. Before first use, carefully inspect the technical condition of the product and make sure that all parts are complete and correctly installed (if applicable, pay particular attention to screws, nuts, rivets and other fastening and structural components). Make sure that the product is not defective or damaged. If any damage or missing parts are found, DO NOT use the product. For further information, contact the seller.

TARGET PATIENT GROUP: Persons suffering from diseases, dysfunctions or injuries for the treatment/rehabilitation or compensation of which this device is intended (see the section concerning the intended purpose of the device in these instructions for use). The product may be purchased independently by the user or on the recommendation of a doctor, therapist or other specialist. Both when purchasing the product independently and when it is recommended by a doctor/therapist/other specialist, the available sizes/required functions and variants of the device, indications and contraindications for use, as well as the information provided by the manufacturer, should be taken into account.

INTENDED PURPOSE: The wrist and thumb orthosis is intended for stabilization and support of the thumb and wrist area. It helps limit excessive movements, improves comfort of use and supports proper positioning of the hand.

INDICATIONS: The wrist and thumb orthosis is indicated for use as support in cases of overuse, thumb pain, wrist discomfort, joint weakness and conditions requiring stabilization of the thumb and wrist. It may be used after injuries, in inflammatory conditions and during activities requiring repetitive hand movements. The product is also suitable for protecting the thumb and wrist during computer work, manual work, physical activity and everyday prevention of overuse.

CONTRAINDICATIONS

Physical or mental limitations (e.g. visual impairment) that prevent safe handling of the product.

The orthosis should not be used in the case of recent injuries, suspected fracture, dislocation, severe sprain or instability of the thumb joint, wrist or hand without prior consultation with a doctor. The product should not be applied directly to open wounds, damaged skin, abrasions, skin inflammation, rash, swelling, infections or other skin lesions at the point of contact with the orthosis. Allergy or hypersensitivity to the materials used in the product is also a contraindication to use. If pain, numbness, tingling, sensory disturbances, impaired circulation, excessive pressure, redness, skin irritation or worsening of symptoms occur, immediately discontinue use of the product and consult a doctor or specialist.

LIST OF COMPONENTS: Wrist and thumb orthosis, instructions for use

PRODUCT FEATURES AT53116: Adjustable hook-and-loop fastening, breathable spacer knit with mesh, lightweight construction

PRODUCT FEATURES AT04303: Adjustable hook-and-loop fastening, removable metal reinforcement, lightweight construction

TECHNICAL SPECIFICATION AT53116

Size: UNIVERSAL

Material: 15% neoprene, 40% nylon, 45% polyester

TECHNICAL SPECIFICATION AT04303

Size: UNIVERSAL

Material: 57% polyamide, 22% polyester, 17% elastane, 4% polypropylene

HOW TO PUT ON THE WRIST AND THUMB ORTHOSIS

Open all hook-and-loop fastenings.

Position the orthosis so that the stabilizing part is on the thumb side and the wide strap surrounds the wrist.

Insert the hand into the orthosis and place the thumb through the designated opening.

Adjust the orthosis to the hand so that it fits closely around the thumb, metacarpus and wrist.

Fasten the strap around the thumb, adjusting the degree of stabilization.

Wrap the wide strap around the wrist and fasten it with the hook-and-loop closure.

Check that the orthosis holds the hand securely but does not cause pain, numbness, tingling or excessive pressure.

If necessary, readjust the fastenings to achieve a comfortable and safe fit.

CAUTION: The orthosis should not restrict circulation or cause discomfort. In the event of pain, numbness, redness or worsening of symptoms, discontinue use.

CAUTION: In the event of a "serious incident" related to the device which directly or indirectly led, might have led or may lead to any of the following events:

- a) the death of a patient, user or other person or
- b) temporary or permanent deterioration in the state of health of a patient, user or other person or
- c) a serious public health threat,

the above-mentioned "serious incident" must be reported to the manufacturer and to the competent authority of the Member State in which the user or patient resides. In Ireland, the competent authority is the Health Products Regulatory Authority (HPRA).

CAUTION: If pain, allergic reactions or other concerning symptoms related to the use of the medical device occur that are unclear to the user, consult a healthcare professional.

CAUTION: It is prohibited to use the product in any manner other than in accordance with its intended purpose.

CLEANING: The orthosis may only be hand washed in lukewarm water (at a temperature of up to 40°C) using a mild detergent (e.g. grey soap). Do not iron. Do not dry clean. Do not use bleach. Protect from moisture. Do not expose to high temperatures. Do not tumble dry. Dry flat, gently squeeze before drying. Note concerning AT04303: Before washing the product, remove the metal reinforcement.

STORAGE: Store the product in a dry place, do not expose it to direct sunlight.

DISPOSAL: After being withdrawn from use, the device should be treated as ordinary municipal waste. Remember to comply with waste segregation regulations. Protect the environment.

SYMBOL EXPLANATION:

Product is flammable – keep away from open flames.

Not suitable for children aged 0–3 years.

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Model:
Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. ANTAR MEDICAL Sp. z o.o. ,03-068 Warsaw, ul. Zawisłańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.
10. The warranty does not cover parts and components subject to natural wear and tear, such as wheels, tires, rubber pads, fuses, etc., the wear of which results from normal use of the product. The warranty period for batteries is 6 months.

ANTAR MEDICAL Sp. z o.o.
03-068 Warsaw, ul. Zawisłańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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